


Sample Name : 120431  
 Acq. Operator : 1  
 Spec. Reported : Integration  
 Method Info :  
 Column: Agilent C18, 100mm\*4.6mm, 3.0um  
 Temperature: 40. Flow rate: 1.5mL/min  
 Mobile phase: Water: ACN 70: 30  
 Tgt Mass (EZX) :  
 Location : Vial 26  
 Inj : 0  
 Inj Volume : 5 ul



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
5.52	0.06	1475.80	528.42	90.08
5.78	0.05	162.65	50.38	9.92