

^1H NMR (400 MHz, $\text{DMSO-}d_6$) δ ppm 1.86 (br. s., 2 H) 2.60 - 2.72 (m, 2 H) 2.76 (dd, $J=13.61, 8.48$ Hz, 1 H) 2.93 (dd, $J=13.68, 5.46$ Hz, 1 H) 3.74 (s, 3 H) 3.98 - 4.13 (m, 1 H) 6.93 - 7.03 (m, 1 H) 7.03 - 7.17 (m, 2 H) 7.20 - 7.40 (m, 1 H) 7.72 (s, 1 H) 7.84 (s, 1 H) 8.37 (br. s., 1 H)

