

^1H NMR (400 MHz, $\text{DMSO-}d_6$) δ ppm 12.75 (s, 1 H) 8.77 (br. s., 3 H) 8.39 (s, 1 H) 7.92 (d, $J=7.58$ Hz, 1 H) 7.44 - 7.53 (m, 4 H) 7.36 - 7.42 (m, 3 H) 7.30 - 7.36 (m, 2 H) 6.96 (d, $J=7.58$ Hz, 1 H) 2.55 (t, $J=7.95$ Hz, 4 H) 2.18 (ddd, $J=11.19, 7.82, 7.64$ Hz, 1 H) 1.78 (ddd, $J=11.43, 8.01, 7.95$ Hz, 1 H)

