


Sample Name : 116050  
 Acq. Operator : 8  
 Spec. Reported : Integration  
 Method Info :  
 Column: Sunfire C8, 150mm\*4.6mm, 3.5um  
 Temperature: 30. Flow rate: 1.0mL/min  
 Mobile phase: Water: ACN 70: 30  
 Tgt Mass (EZX) :  
 Location : P2-F-06  
 Inj : 0  
 Inj Volume : 2.6 ul



Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
7.47	0.13	17.65	4.60	0.22
7.92	0.15	7879.18	1774.37	98.29
8.40	0.13	119.36	27.98	1.49