



Sample Name : 146596
Acq. Operator : C
Spec. Reported : MS Integration
Tgt Mass (EZX) :
Location : P1-D-07
Inj : 1
Inj Volume : 6 ul




Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
2.47	0.06	9391.38	2606.12	97.94
2.63	0.05	197.47	68.71	2.06

Ret. Time: 2.53

<<<< POSITIVE SPECTRA >>>>

ES-API Positive



Ret. Time: 2.69

ES-API Positive

