

$^1\text{H}$  NMR (400 MHz, CHLOROFORM-*d*)  $\delta$  ppm 6.94 - 7.02 (m, 1 H) 7.02 - 7.09 (m, 1 H) 7.19 - 7.24 (m, 2 H) 7.29 - 7.36 (m, 1 H) 7.43 - 7.51 (m, 2 H) 7.51 - 7.60 (m, 1 H) 8.08 (dd,  $J=8.09, 1.67$  Hz, 1 H) 10.50 (s, 1 H)

