

^1H NMR (400 MHz, $\text{DMSO-}d_6$) δ ppm 1.39 (d, $J=10.15$ Hz, 2 H) 1.45 - 1.59 (m, 2 H) 1.94 (t, $J=10.82$ Hz, 2 H) 2.72 (d, $J=11.00$ Hz, 2 H) 3.48 (s, 2 H) 3.72 - 3.89 (m, 1 H) 4.39 (s, 2 H) 7.11 - 7.29 (m, 4 H) 7.29 - 7.43 (m, 6 H) 7.47 (d, $J=7.46$ Hz, 1 H) 7.53 (t, $J=7.70$ Hz, 1 H) 7.61 (d, $J=7.95$ Hz, 1 H) 7.66 - 7.77 (m, 2 H) 7.83 (s, 1 H) 8.34 (s, 1 H) 8.37 (dd, $J=4.10, 2.14$ Hz, 1 H)

