



Sample Name : 147920
 Acq. Operator : C
 Spec. Reported : MS Integration

Seq. Line : 0
 Location : P1-E-02
 Inj : 0
 Inj Volume : 8 ul

MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



MWD1 A, Sig=214,4 Ref=off



Integration Results for MWD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%	MS(+)
2.25	0.06	158.18	44.31	0.60	*** 38
2.55	0.21	25815.21	1774.10	98.37	*** 50
3.42	0.04	270.55	103.32	1.03	*** 50

Ret. Time: 2.60

<<<< POSITIVE SPECTRA >>>>

