


Sample Name : 114739 Tgt Mass (EZX) :
 Acq. Operator : 1 Location : P1-F-01
 Spec. Reported : Integration Inj : 0
 Method Info : Inj Volume : 5 ul

Column: Agilent C18, 100mm*4.6mm, 3.0um
 Temperature: 40. Flow rate: 1.5mL/min
 Mobile phase: Water: ACN 70: 30



Integration Results for DAD1 A, Sig=210,4 Ref=off

RetTim	Width	Area	Height	Area%
3.35	0.18	3425.10	583.62	99.84
4.28	0.05	5.42	1.81	0.16