


Sample Name : 131365
 Acq. Operator : 1
 Spec. Reported : Integration
 Method Info :
 Column: Sunfire C8, 150mm*4.6mm, 3.5um
 Temperature: 40. Flow rate: 1.5mL/min
 Mobile phase: Water: ACN 70: 30

Tgt Mass (EZX) :
 Location : P1-E-06
 Inj : 0
 Inj Volume : 6 ul



Integration Results for DAD1 A, Sig=223,4 Ref=off

RetTim	Width	Area	Height	Area%
7.63	0.12	3697.68	1014.00	99.92
9.16	0.05	3.11	0.79	0.08