

¹H NMR (500 MHz, DMSO-*d*₆) δ ppm 9.06 (s, 1 H) 7.76 (d, *J*=1.96 Hz, 1 H) 7.64 - 7.71 (m, 2 H) 7.60 (dd, *J*=8.56, 2.20 Hz, 1 H) 7.54 (d, *J*=8.31 Hz, 1 H) 7.42 - 7.51 (m, 2 H) 7.36 - 7.42 (m, 1 H) 5.15 (dd, *J*=8.80, 3.42 Hz, 1 H) 4.64 (d, *J*=10.76 Hz, 1 H) 4.46 - 4.55 (m, 2 H) 4.38 - 4.46 (m, 1 H)

