

^1H NMR (400 MHz, $\text{DMSO-}d_6$) δ ppm 2.37 (s, 3 H) 2.43 (s, 3 H) 2.71 - 2.80 (m, 2 H) 2.82 - 2.96 (m, 2 H) 5.19 (s, 2 H) 7.00 (m, $J=8.11$ Hz, 2 H) 7.19 (m, $J=8.11$ Hz, 2 H) 7.50 - 7.59 (m, 2 H) 7.59 - 7.76 (m, 2 H)

