

¹H NMR (400 MHz, DMSO-*d*₆) δ ppm 3.20 - 3.26 (m, 1 H) 3.30 (br. s., 1 H) 3.43 - 3.54 (m, 3 H) 3.71 (dd, *J*=10.58, 5.20 Hz, 1 H) 4.62 (d, *J*=5.62 Hz, 3 H) 5.09 (d, *J*=5.26 Hz, 1 H) 5.13 (dd, *J*=4.83, 2.26 Hz, 2 H) 5.17 (d, *J*=7.70 Hz, 1 H) 5.56 (t, *J*=5.81 Hz, 1 H) 7.28 (s, 1 H) 7.66 (d, *J*=1.47 Hz, 1 H) 7.71 (dd, *J*=7.95, 1.71 Hz, 1 H) 7.82 - 7.99 (m, 2 H) 12.88 (s, 1 H)

