



Sample Name	: 146366	Location	: P1-C-01	Seq. Line	: 0
Acq. Operator	: C			Inj	: 0
Spec. Reported	: MS Integration			Inj Volume	: 7 ul

MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



MWD1 A, Sig=214,4 Ref=off



Integration Results for MWD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%	MS(+)
2.68	0.09	15127.31	2547.93	100.00	*** 57

Ret. Time: 2.73

<<<< POSITIVE SPECTRA >>>>

